

# The Warwick Range


The Warwick Range from Warwick SASCo is a wide range of sterilizable long-life plastic polyware products.

Designed as an alternative to stainless steel products, they are ideal for use in Decontamination Services, Operating Rooms, Clinics, Surgeries and other surgical and healthcare services.



# The Warwick Range

Warwick products are manufactured from a high quality polypropylene material that is approved for use in medical devices. All products in the Warwick Range are registered as CE Class 1 Medical Devices.



Kidney Dishes and Gallipots • Lotion Bowls and Washbowls • Instrument Trays  
Instrument Protection and Transportation • Scope and Instrument Disinfection Sets  
Silicone Protection Mats • Traceability Tags • Measuring Jugs

## Warwick Benefits

- Low cost.
- Robust and hard wearing - guaranteed for 1000 washes.
- Polypropylene is easier to clean than stainless steel as blood protein is easily removed from the surface, leaving no residue.
- 50% lighter than stainless steel, assisting users with storage, lifting and handling.
- Can be sterilized using Autoclave, ETO or Hydrogen Plasma.
- Can be used as stand-alone items or as a support component for surgical sets to assist surgical and hospital teams with surgical procedures.
- Widely used in Sterile Services for packing instrument sets and procedure kits for the Autoclave.
- Disinfection sets can be used for manual sterilization in Sterile Services, Clinics, etc.
- Can also be used non-sterile in hospital wards and clinics as vomit bowls, collection bowls.
- All Warwick products have smooth corners and surfaces, reducing risk of paper rips and the need for re-processing.

**WarwickSASCO**  
Moulding The Future of Medical Plastics